



Liberty™ 2.0 Series

Automated Microwave Peptide Synthesizers

Introducing Our Liberty 2.0 Series


The Liberty Blue™ 2.0 and Liberty PRIME™ 2.0 are the most efficient automated microwave peptide synthesizers available. Released in 2022, this new generation of systems extends beyond the award-winning High-Efficiency SPPS (HE-SPPS) methodology and instrumentation introduced by CEM's best-selling Liberty Blue synthesizer. The Liberty 2.0 Series of instruments features significant advancements in peptide quality, robustness, and overall flexibility compared to the previous generation. The biggest breakthrough offered by these instruments is the ability to substantially decrease solvent use while increasing the complexity of peptides that can be synthesized with high purity. Engineering advancements to system fluidics, reaction vessel capabilities, including Headspace Flushing, and software provide revolutionary robustness and flexibility. Learn more about the exciting updates to the most advanced peptide synthesizer line available.



Improved Peptide Quality

The Liberty Blue 2.0 and Liberty PRIME 2.0 provide the ability to make even higher purity peptides than previously possible with their respective predecessors. This is based on advances in the system fluidics, reaction vessel environment, and chemistry methodology.

- Optimized chemistry methodology with ProTide resins.
- Cleanest reaction environment with updated headspace flushing technology.
- More direct optimized flow paths.
- Digitally variable N₂ control for optimized bubbling, purging, and draining.



Increased Reliability



All Liberty 2.0 level systems feature more direct and cleaner flow paths. The result is a new level of robustness providing higher runtime with less downtime even under high usage conditions.



Ultimate Flexibility

You can do almost anything with the Liberty 2.0 series.


- Access a wide scale range from 5 μ mol – 5 mmol.
- Wider pressure range for accessing viscous green solvents (TamiSolve).
- Automate synthesis of complex steps and more.



Expanded Molecular Access

Protein synthesis on the Liberty PRIME 2.0 demonstrates the enhanced performance of the 2.0 series to synthesize very long peptides and even small proteins. The combination of the chemical methodology developed by CEM with the engineering improvements of the new 2.0 series push the limits of chemical synthesis.


Ubiquitin



Length: 76 AA

Synthesis Time: 9 hrs. 35 min.


HIV Protease



Length: 99 AA

Synthesis Time: 12 hrs. 48 min.

MDM2



Length: 118 AA

Synthesis Time: 16 hrs. 26 min.

Collagen Like Peptide




Length: 99 AA

Synthesis Time: 12 hrs. 30 min.

Headspace Flushing — Cleaner Reaction Vessel Surfaces


A new feature with the Liberty Blue 2.0 and Liberty PRIME 2.0 systems is the use of headspace flushing to provide a cleaner reaction vessel surface from volatile deprotection bases at elevated temperature. The use of headspace flushing is beneficial for achieving unprecedented levels of purity for longer sequences.



Note: Headspace flushing technology is utilized on the Liberty Blue 2.0 and Liberty PRIME 2.0 systems

Synthesize Peptides with High Efficiency

The Liberty Blue 2.0 is the exciting next generation of the best-selling peptide synthesizer in the world. Compared to its predecessor, the Liberty Blue 2.0 features the ability to make cleaner and longer peptides with a new level of robustness. The system features a 4-minute cycle time along with a 90% solvent reduction based on High Efficiency Solid Phase Peptide synthesis (HE-SPPS). The Liberty Blue 2.0 features the latest in engineering for fluidic deliveries, true internal temperature feedback control, and software control with 21 CFR Part 11 compliance.






Unparalleled Speed and Efficiency

The Liberty PRIME 2.0 is the absolute best peptide synthesizer available in the world. The system features a unique one-pot deprotection and coupling process allowing for a remarkable 2.5 minute cycle time with only 8.5 mL waste per cycle (at 0.1mmol). A key improvement in the 2.0 version of the Liberty PRIME is the ability to make cleaner and longer peptides with a new level of robustness than previously possible.

The Liberty PRIME 2.0 provides the ultimate flexibility for long, difficult, or complex peptides. Additionally, the system provides an ultra-fast sequential synthesis workflow that can dramatically outperform traditional parallel peptide synthesizers. The Liberty PRIME 2.0 can be utilized in cGMP peptide production.

- Massive reduction in solvent and deprotection usage.
- Dramatically faster deprotection step because the reaction mixture remains at high temperature.
- Shorter overall cycle time due to the removal of the post-coupling draining step.



Find Your Automated Needs

All Liberty synthesizer models give you access to:

- Fast, Flexible, High-Purity Synthesis
- Synthetic Biomolecules
- Alternative Solvents
- Complex side-chain functionalization, cyclization, and branching
- High-throughput synthesis
- Less Waste

	Liberty Blue 2.0	Liberty PRIME 2.0
		
Cycle Time (at 0.1mmol)	4 minutes	2.5 minutes
System Waste (at 0.1mmol)	16 mL	10 mL
Scale Range	0.005 – 5 mmol	0.005 – 5 mmol
Headspace Flushing	Yes	Yes
Amino Acid Positions	27	30
External Reagent Positions	4	5
RV Camera	Yes	Yes
LED Visual Feedback	Yes	Yes
High-Throughput Options	HT4, HT12	HT4, HT12, HT24
Optional Accessories	N/A	Full cGMP Compliance Package




Liberty PRO









Easily scale up peptides optimized on any Liberty 2.0 series system. The Liberty PRO™ allows for fully automated, production scale, microwave peptide synthesis under cGMP using a 3, 8, or 15 liter reaction vessels. Batch crude peptides up to 1 kg can be synthesized with a typical cycle time of 15 – 45 min per amino acid.



No Washing — Solid-Phase Peptide Synthesis (UE-SPPS)

Traditional SPPS generates large amounts of waste from repetitive washing after every coupling and deprotection step (about 80% of total waste). The Liberty Blue 2.0 and Liberty PRIME 2.0 systems eliminate all washing during each cycle! This is based on improved chemical methodology and a new headspace flushing technology for removing volatile deprotection base at high temperatures. Run High Efficiency SPPS (HE-SPPS) with limited waste or UE-SPPS with absolutely no wash waste on the Liberty 2.0 series of automated microwave synthesizers to access standard, complex, and very long peptides at high purity and yield.






Unparalleled Peptide Synthesis




The Most Powerful Partnering of Chemistry and Engineering

The Liberty 2.0 series, like preceding models, utilizes CEM's High-Efficiency Solid Phase Peptide Synthesis (HE-SPPS) and CarboMAX coupling methodologies.




Pioneers in Microwave SPPS

- Rapid, responsive heating
- Energy efficient
- Optimized methodology for high speed and purity



One-pot Coupling/Deprotection


CEM's recently developed one-pot coupling/deprotection methodology provides unparalleled simplicity and efficiency, specifically featured on the Liberty Prime 2.0.





Automated High Throughput Production

The Liberty series of instruments enable high throughput peptide synthesis with high purity. Optional High-Throughput modules allow for running automated batch sizes of 4 (HT4), 12 (HT12), or 24 (HT24) peptides on a Liberty 2.0 series system.



Sequence	Crude Purity	Final Purity*	Recovery (%)	mg
GWVKPIIGHHAYGDQYRAT	73	99.1	61	103.2
TLYEQEIEV	49	94.5	68	39.7
HGSRKNITDMVEGAKKANG	73	91.7	57	91.1
SLLNQPKAV	79	99.7	46	35.0
EDPYLFELPVLKYLDMGTT	76	98.5	68	38.4
ALAVLSNYDA	84	99.3	13	11.4
TMEDKIYDQQVTKQCLCF	47	97.1	44	47.2
YSYPETPLYMQTASTSYE	47	95.2	26	29.7
KVGYTERQRWDFLSEASIM	61	97.3	44	63.2
RLRMREHMMKNVDNQNQD	65	96.3	51	75.2
VYEKNGYIYF	90	98.5	47	54.8
ALAVLCNYDA	73	99.3	45	34.7
ALVPPSKRKMWVVSPEAKA	78	97.1	66	111.2
ISTPTPTIVHPGSLPLHLG	75	99.2	56	81.7
IVQENNTPGTYLLSVSARD	74	96.3	56	89.7
RFHMKVSVYLLAPLREALS	75	96.8	51	88.7
ENLKQNDISAEFTYQTKDA	82	99.3	61	111.5
YMPVNSEV	70	97.8	44	33.3
TNDVKTLDLNGVIEEFT	59	97.6	24	29.9
SAWLFRMWYIFDHNYLKPL	48	99.9	66	79

*The peptides were purified on the Prodigy to give >90% final purity on all samples.


Hilf, N. et al (2019) Nature 565, 240-245.



Flexible and Powerful Software


The Liberty 2.0 Series of Peptide Synthesizers feature an updated software system that provides the latest in ease-of-use, flexibility, and conformity with 21 CFR Part 11 compliance. The software allows for simple customization of reagent positions, peptide modifications, and the ability to edit/stop during a method run.

Easy and Flexible Sequence Programming




The method editor on the Liberty 2.0 series allows full customization of reaction parameters. Create, edit and store methods to streamline the run process.

Program Post-Synthesis Modification




Versatile options for post-synthetic steps like branching, cyclization, orthogonal deprotection, labeling, and more. Empower your chemistry with flexible options.

Full Suite of Calculator Tools




Convenient calculators provide the reagents and amounts for a run, streamlining the process of reaction preparation.

Edit or Monitor at any Step



Easily edit and monitor reactions from the run screen. With an HT module, program a queued run, even while a synthesis is in progress.



Liberty Methods

Name: Liraglutide

Tags: +

Sequence (N → C)

A Ala 1	R Arg 2	N Asn 3	D Asp 4	C Cys 5	D Glu 6	E Glu 7	D Gly 8	H His 9	I Ile 10
L Leu 11	K Lys 12	M Met 13	F Phe 14	P Pro 15	S Ser 16	T Thr 17	W Trp 18	Y Tyr 19	V Val 20
1 EX1 21 (EX1)	2 EX2 22 (EX2)	3 EX3 23 (EX3)	4 EX4 24 (EX4)	5 EX5 25 (EX5)	6 EX6 26 (EX6)	7 EX7 27 (EX7)	8 EX8 28 (EX8)	9 EX9 29 (EX9)	10 EX10 30 (EX10)

Resin Name:

Resin Loading: 0

Final Deprotection Cycle: 0 To Final Deprotection

Resin Cycle: 0 To Resin Swelling

Index	Amino Acid	Cycle
-------	------------	-------

Forget changes

Save

Paste sequence

First Synthesis Step

Load default cycles

Components

Fluidics	Connected
Pressure	15.00 PSI
Microwave	Connected
Temperature	50 °C
PRIME Pump	Connected


Sensors

Waste	OK
SpS	OK
Volume SpS	OK
RF Overflow	OK

Maintenance

Reagent	OK
Wash	OK
50%	OK
50%	OK
12%	OK

We have You Covered from Spot Synthesis to Kilogram Scale Production



Cleavage

Razor

The RAZOR® reduces standard cleavage times from 3 – 4 hours down to only 30 min. The system performs cleavage of up to 12 peptide resins at a time with precise temperature control. The system is compact and easily fits in standard fume hoods.



Purification

Prodigy

The Prodigy™ Peptide Purification System is a Preparatory HPLC that is specially engineered for peptide purification. It features elevated temperature purification, a focused gradient calculator, and is capable of purifying milligrams to grams of crude peptide — all with modern and easy to use software.





We Simplify Science

cem.com



United States (Headquarters): 800-726-3331 | info@cem.com
For distributors and subsidiaries in other regions, visit cem.com/contact